강음이온 교환 멤브레인 Adsorber(크로마토그래피)
(Purcise™ Q Strong Anion Exchange Membrane Adsorber (Chromatography))

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멤브레인 사양(Membrane Details)

Product Series	**Purcise™ Q**
Membrane Materials	PES
Functional Group	Mixed amino group
Typical Operating Flow Rate	1-25 MV/min
Nominal Pore Size	0.7 μm
BSA Binding Capacity	>50 mg/mL
DNA Binding Capacity	>20 mg/mL
Pre-use Storage Conditions	Store at room temperature or 2~8°C
Post-use Storage Solution	0.1M NaOH

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Specification	CXU33	CXD32	
Lab Scale			
Volume	0.2 mL	0.45 mL	0.9 mL
Minimum standards BSA binding capacity	50 mg/ml	50 mg/ml	50 mg/ml
Membrane Configuration	Flat sheet	Flat sheet	Flat sheet
Layers	3	8	16
Connection forms	I : Female Luer; O: Similar to Male Luer	Female Luer	Female Luer
Structural material	Upper : PP Lower : PP	Upper : PP Lower : PP	Upper : PP Lower : PP
Recommended flow rate range*	0.2-5 mL/min	0.45-11.25 mL/min	0.9-22.5 mL/min
Max. operating pressure	4.0 bar	4.0 bar	4.0 bar
Typical application	Process development experiments; Scaled down laboratory tests		

Specification	CX0005	CX0010	CX0025	CX0050
Lab Scale				
Volume	5.0 mL	10.0 mL	25.0 mL	50.0 mL
Minimum standards BSA binding capacity	50 mg/ml	50 mg/ml	50 mg/ml	50 mg/ml
Membrane Configuration	Pleated	Pleated	Pleated	Pleated
Layers	8	16	8	16
Connection forms	M6 Tapered Thread Port	M6 Tapered Thread Port	TC25	TC25
Structural material	Membrane Support Layer / Drainage Structure: PP Core / Support / End Cap: PP Membrane Shell: PP Seal Ring: Red Silicon Valve: PP			
Recommended flow rate range*	5-125 mL/min	10-250 mL/min	25-625 mL/min	50-1250 mL/min
Max. operating pressure	3.0 bar	3.0 bar	3.0 bar	3.0 bar
Typical application	Process Development Experiment	Process Development Experiment	Pilot-Scale Experiment	Pilot-Scale Experiment

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
Specification	CX0075	CX0150	CX0225	CX0450
Pilot & Production Specification				
Volume	75 mL	150 mL	225 mL	450 mL
Minimum standards BSA binding capacity	50 mg/ml	50 mg/ml	50 mg/ml	50 mg/ml
Membrane Configuration	Pleated	Pleated	Pleated	Pleated
Layers	8	16	8	16
Connection forms	TC25	TC25	TC25	TC25
Structural material	Membrane Support Layer /Drainage Structure: PP Core / Support / End Cap: PP		Membrane Shell: PP Seal Ring: Red Silicon	Valve: PP
Recommended flow rate range*	0.075-1.875 L/min	0.15-3.75 L/min	0.225-5.625 L/min	0.45-11.25 L/min
Max. operating pressure	3.0 bar			

Specification	CX1000	CX2000	CX2500	CX5000
Pilot & Production Specification				
Volume	1000 mL	2000 mL	2500 mL	5000 mL
Minimum standards BSA binding capacity	50 mg/ml	50 mg/ml	50 mg/ml	50 mg/ml
Membrane Configuration	Pleated	Pleated	Pleated	Pleated
Layers	8	16	8	16
Connection forms	TC50	TC50	TC50	TC50
Structural material	Membrane Support Layer /Drainage Structure: PP Core / Support / End Cap: PP		Membrane Shell: PP Seal Ring: Red Silicon	Valve: PP
Recommended flow rate range*	1-25 L/min	2-30 L/min	2.5-62.5 L/min	5-75 L/min
Max. operating pressure	3.0 bar			

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
Laboratory scale-up performance
(in terms of BSA-binding capacity)

Actual Test Data for Purcise™ Q Removal of Host Protein (HCP) from Antibody Samples

Actual Test Data 1				
Samples	HCP	HMW	SEC Main Peak	LMW
Initial samples	28 ng/mg	2.07 %	97.89 %	0 %
pH 8.2 Cond.~4 mS/cm	3.5 ng/mg	1.73 %	98.26 %	0 %
pH 7.75 Cond.~4 mS/cm	4.8 ng/mg	1.65 %	98.35 %	0 %
pH 7.5 Cond.~4 mS/cm	4.3 ng/mg	1.59 %	98.41 %	0 %

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Actual Test Data 3			
Samples	Description	Load	HCP
Initial samples	Initial samples	N/A	216 ppm
Samples 1	pH 7.5, Cond. <5mS/cm through hybrid samples	0.5 g/ml	< 9 ppm
Samples 2	pH 7.5, Cond. <5mS/cm through hybrid samples	2.0 g/ml	15 ppm
Samples 3	pH 7.5, Cond. <5mS/cm through hybrid samples	4.0 g/ml	20 ppm
Samples 4	pH 7.5, Cond. <5mS/cm through hybrid samples	4.0 g/ml	24 ppm

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Challenge Test Conditions

Parameter	Condition 1 - Two Parallel Groups	Condition 2 - Two Parallel Groups
pH and Conductivity	pH 7.0, Conductivity ~6 mS/cm	pH 8.0, Conductivity ~6 mS/cm
Membrane Volume	0.9 ml (Lot. CXD32EAQ16CC1P)	
Flow Rate	10 MV/min 9 ml/min	
Challenged Virus Types	MVM, X-MuLV, Reo-3, PrV	
mAb Concentration	9.42 g/L	
Load Volume	3.1 kg/L (g/ml)	
Residual HCP in Flow Through	~20 ppm	

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엔도톡신 제거 테스트(Purcise™ Q Endotoxin Removal Test)

Test Conditions

Conducting DOE Experiments to Study the Performance of Q Membrane in Removing Endotoxins from Formulation Buffers under Low Load Conditions

Condition	Low	High
pH	4	9
Sucrose Concentration	0 % (W/V)	15 % (W/V)
NaCl Concentration	0 mM	200 mM

Test solution spiked with endotoxin concentration: 100 EU/ml
Purcise™ Q membrane volume: 0.9 ml (configured with 16 layers)
Purcise™ Q membrane filtration spiked solution volume: 90 ml
Filtration flow rate: 9 ml/min (10 MV/min)
Collect the solution before and after Q membrane filtration, and measure the endotoxin concentration. The endotoxin removal is expressed using log reduction value (LRV).

Test Results


Endotoxin Clearance LRV =1.60 + 0.2675*A - 0.5826*B - 0.6498*C + 0.2937*BC + 0.8019*A^2 + 0.3434*B^2 + 0.3048*C^2, where

A: pH

B: Sucrose Concentration

C: NaCl Concentration

Response Surface Curves

At **pH 4**, the response surface curve of B: sucrose concentration and C: sodium chloride concentration.

At **pH 6.5**, the response surface curve of B: sucrose concentration and C: sodium chloride concentration.

At **pH 9**, the response surface curve of B: sucrose concentration and C: sodium chloride concentration.

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Copyright © BVSKOREA.CO., LTD. All Rights Reserved