
Consieve® UET

Consieve® UET Cassette User Manual

Beyond Valuable Solution

Please contact us for more information

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Copyright @ BVSKOREA.CO.LTD. All Rights Reserved

1. Cassettes

Cassette boxes include the following components:

- Cassette in sealed plastic bag
- Cassette gasket
- Certificate of compliance
- User manual

Cassette structure

❶ Inlet (feed) port
❷ Outlet (retentate)
❸ port Permeate port
❹ Cassette information
❺ Alignment notches
❻ Incorporated gasket


Figure 1. The main parts of cassette

Product labeling

The cassette label is laser-engraved on the side of the cassette.
It contains manufacturer, lot number, cassette type (Lab or Flow), part number, nominal molecular weight cutoff (NMWC), and membrane area in square meters and square feet (Fig 2).

Figure 2. Sample product label

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1. Remove the tie rod nuts, washers,and spacers. Slide the backing plate away from the flow distribution manifold (Figure 3).
2. To install one cassette, clean and wet a silicone gasket by rinsing it with DI water or water-for-injection (WFI). Place the gasket against the flow distribution manifold, aligning the holes in the gasket with the holes in the manifold. The gasket will stick to the manifold. Place the cassette into the holder with the gasket-side of the cassette facing the backing plate, and slide the cassette against the gasket. To install two or more cassettes, complete the step above for installing one cassette.Place the second cassette into the holder between the first cassette and the backing plate. Ensure the gasketed side of the second cassette faces the backing plate. Place subsequent cassettes in the holder in similar fashion.
3. Slide the backing plate forward to hold the cassettes in the hardware, ensuring the holes in the manifold, gasket, and cassettes line up.
4. Add spacers to the tie rods until about 0.75-inches (16mm) of thread is exposed on each tie rod. Reinstall the washers and nuts. Finger tightened the nuts evenly.


Figure 3. Cassette Holder

❶ Flow distribution manifold
❷ Backing plate
❸ Gasket
❹ Cassette
❺ Tie rod
❻ Spacer
❼ Washer
❽ Nut
❾ Feed
❿ Retentate
⓫ Permeate

5. Tighten each nut alternately 1/4 turn at a time with the supplied torque wrench and socket until each nut is tightened to the values listed in Table 1 above. If a seal is not achieved, there may be a cassette, gasket, or surface finish irregularity. Discontinue torque tightening and inspect the sealing surfaces.

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3. Using a cassette

Introduction

This section provides basic procedures for cassette preparation, operation, cleaning, and storage. This section does not cover process optimization.

Cassette preparation

Before using a new or used-and-stored cassette, rinse the storage solution from the cassette following these steps:

1. Open the cassette bag, drain the cassette, and rinse the storage solution from the exterior of the cassette with high purity water. Use appropriate personal protective equipment to avoid contact with the NaOH/glycerin storage solution.
2. Configure a filtration system according to Figure 4. Install the cassette in the cassette holder. (See the installing a cassette for specific installation instructions.) Close the feed, retentate, and permeate valves.
3. Direct the retentate and permeate lines to waste.
4. Fill the reservoir with deionized (DI) water or water for injection (WFI) (5 liters per 0.09 m^2 [1.0 ft^2] of membrane area).
5. Open the feed and retentate valves. Close the permeate valve. Pump about 10% of the water through the retentate line to waste at a feed pressure of 0.3 barg (5 psig).


Figure 4. Typical setup for rinsing cassettes

6. Open the permeate valve and close the retentate valve. Pump the remaining water through the permeate line to waste at a transmembrane pressure (TMP) of 0.3 barg (5 psig).
7. Optionally, analyze a permeate sample to verify sufficient flushing for your specific application. Flush with additional water, if needed.

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Measuring water flux

Measuring water flux involves measuring the flow of water through the cassette's membrane under controlled conditions. You convert the measurement into units that represent water flux. The water flux value you obtain provides an indicator of the performance of the cassette.

1. Setup your cross-flow system to circulate the retentate and permeate back to the feed reservoir (Figure 5).
2. Fill the feed reservoir with a volume of deionized water or WFI. Fully open the feed, retentate, and permeate valves.
3. Start the feed pump and adjust it and the retentate valve to obtain a transmembrane pressure of 0.3 barg (5 psig). Using a graduated cylinder and stopwatch, measure and record the permeate flow rate in mL/min. Record the temperature of the water in the feed reservoir.
4. Adjust the feed pump and retentate valve to produce a transmembrane pressure of 1 barg (15 psig). Using a graduated cylinder and stopwatch, measure and record the permeate flow rate in mL/min. Record the temperature of the water in the feed reservoir.


Figure 5. Typical configuration for measuring water flux

5. Drain the system.
6. Calculate flux. Convert the permeate flow rate measurements to flux using this equation:

Flux in LMH (l/m²/hr) = {permeate flow in ml/min ÷ cassette surface area in square meters} × 0.06

7. Plot the calculated flux values against the transmembrane pressure (Figure 6). Determine the water flux in LMH at 0.7 barg (10 psig). In Figure 6, the flux at 0.7 barg (10 psig) is 115 LMH.


Figure 6. Plot of water flux versus transmembrane pressure

8. Multiply the calculated flux value by the appropriate viscosity correction factor in Table 2. For example, if the flux at 0.7 barg (10 psig) is 115 LMH, and the temperature of the water during permeate flow measurement was 18 °C, then the flux normalized to 20 °C is:

115 LMH x 1.053 = 121 LMH

9. To standardize the flux value to 1 psig, divide the normalized water flux value by 10 psig. For example,

121 LMH ÷ 10 = 12.1 LMH/psig

10. Record the normalized flux value in your laboratory notebook or production records as a benchmark for future comparison.

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Table 2. Viscosity correction factor

Temperature in °C when permeate flow was measured	Viscosity correction factor	Temperature in °C when permeate flow was measured	Viscosity correction factor
T (°C)	F	T (°C)	F
4	1.567	25	0.89
5	1.519	26	0.871
6	1.472	27	0.851
7	1.428	28	0.833
8	1.386	29	0.815
9	1.346	30	0.798
10	1.307	31	0.781
11	1.271	32	0.765
12	1.235	33	0.749
13	1.202	34	0.734
14	1.169	35	0.719
15	1.139	36	0.705
16	1.109	37	0.692
17	1.081	38	0.678
18	1.053	39	0.665
19	1.027	40	0.653
20	1	41	0.641
21	0.978	42	0.629
22	0.955	43	0.618
23	0.933	44	0.607
24	0.911	45	0.596

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Integrity test

Before using your cross-flow system and cassette to process product, test the cross-flow system to ensure itdoes not leak. Then test the integrity of the cassette.
To test your cross-flow system for leaks, pressurize it with water and check the tubing or piping connections for leaks. To test the integrity of the cassette, pressurize the system with compressed air and measure the diffusionalflow of air through the cassette membrane. A high flow rate (above the cassette's specification) indicates there might be a hole or tear in the membrane.

1. Ensure the compressed air valve is closed and the pressure regulator is set to 0 barg (0 psig). Close the feed and drain valves, and open the retentate and permeate valves.
2. Remove water from the cassette feed-retentate path by opening the compressed air valve, and slowly adjusting the pressure regulator to 0.3 barg (5 psig). Flow air through the system until water stops draining from the retentate line. Close the retentate valve to allow the air pressure to drive the water from the permeate port and permeate line.
3. Slowly adjust the pressure regulator to 1 barg (15 psig).
4. Fill a raduated cylinder with water and place it inverted in a beaker half filled with water.Connect flexible tubing to the permeate outlet (Figure 7).
5. When the air bubble rate is stable, note the time and the standing air volume in the graduated cylinder.


Figure 7. Configuration for integrity test

6. Using the time and air volume data, calculate the air diffusion rate in ml/min@1bar(15psi)
7. Compare your results with the specifications below. If your measured flow rate is less than the specification, the cassette is acceptable for use.

Table 3. Integrity Test Specification

Lab Packet	Lab	Flow-5	Flow-25
≤ 1.1ml/min	≤ 12ml/min	≤ 50ml/min	≤ 250ml/min

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Lab Packet	Lab	Flow-5	Flow-25
≤ 1.1ml/min	≤ 12ml/min	≤ 50ml/min	≤ 250ml/min

Operating conditions

Drain the water. Add a sufficient volume of buffer to the feed tank to maintain an adequate working volume. Circulate buffer through the filtration system for up to five minutes to condition it prior to addition of process fluid. Drain the buffer.

When using cassettes, process optimization is highly recommended. Typically, an effective cross flow rate will range from 5.4 to 8.6 liters per m2 (500-800 ml per ft^2) of membrane, and the permeate pressure will be zero. The maximum feed pressure for cassettes is 4 barg (60 psig) at 50°C.

Before beginning a filtration, it may be useful to generate a pump curve that relates pump speed to flow rate. This curve allows easy adjustment of filtration flow rate, eliminating the need to physically measure the flow rate and adjust through trial and error. Alternatively, a flowmeter may be used.

Product recovery

The yield of product in the retentate is generally increased by using a product recirculation after filtration.

1. Close the permeate valve.
2. Reduce the cross flow to 10% of the process cross flow.
3. Recirculate the fluid for 5-15 minutes.

To recover product from the feed tank, use the drain valve that, depending on system configuration, may be before or after the pump. To evacuate product from the remainder of the system, specific system design will dictate whether to operate the pump slowly in either the forward or reverse direction or to drain the system by gravity.

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Cassette cleaning

Cassettes should be cleaned properly before reuse or storage. Before cleaning, cassettes should be rinsed with buffer and water.
It is highly recommended to optimize the cassette cleaning methods. The following is a general cassette cleaning protocol:

1. Direct the retentate and permeate lines to the feed tank (Figure 8). Open the feed and retentate valves. Close the permeate valve.
2. Circulate a minimal volume of buffer across the retentate side for 5-10 minutes at 1-1.5× the process cross flow.
3. Drain the system. Add a minimal volume of water and circulate water across the retentate side for 5-10 minutes at 1-1.5× the process cross flow.
4. Drain the system. Add cleaning solution to the feed tank at a ratio of 15-20 liters per m^2 (1.5-2 liters per ft^2) of membrane area. See Table 4 for typical cleaning conditions.
5. Open the feed and retentate valves and pump about 10%of the cleaning solution through the retentate line to waste.
6. Open the permeate valve, close the retentate valve, and pump about 10% of the cleaning solution through permeate line to waste.
7. Stop the pump and direct the retentate and permeate lines into the feed tank. Open the feed, retentate, and permeate valves.
8. Start the pump and adjust the cross flow to 1-1.5× the process cross flow. Circulate the cleaning solution for 30-60 minutes.
9. Drain the cleaning solution from the system.
10. Flush the system using 22-55 liters of clean water per m^2 (2-5 liters per 1.0 ft^2) of membrane area and the same process described in steps 5 through 9 above. The circulation time may be reduced to 15 minutes.

Cleaning efficacy is commonly evaluated by comparing clean water fluxes pre- and post-use. It is good practice to record clean water flux after each cleaning. Commonly, users will discard cassettes if the water flux has reduced to 60-80% of original water flux. Note that reduced water flux often does not indicate that process fluid flux will likewise be reduced.


Figure 8. Typical setup for recirculation of cleaning solution

Table4. Recommended cleaning conditions

Cleaning agent	Cleaning conditions Contact
0.1-0.5 N sodium hydroxide	time = 60 minutes Temperature = 40°C (113°F)
0.5 N NaOH with 100-300 ppm sodium hypochlorite	Contact time = 60 minutes Temperature = 20°C (68°F)

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Storing a cassette

Cleaned cassettes should be stored according to the conditions outlined in Table 5. These conditions keep the membrane wetted and inhibit microbiological growth. Introduce the storage solution into the cassette as per the instructions for cleaning, steps 4 through 8. The recirculation time (step 8) may be reduced to 10-15 minutes.

It is good practice to store cassettes under refrigeration, and this is especially recommended for long-term storage (i.e., 1-12 months). Store cassettes in air-tight containers and add a small volume of storage solution to the containers to minimize evaporation. Alternatively, cassettes may be submerged in storage solution for up to one month. For short-term storage (up to 3 days), cassettes may be stored in holders that are compressed through a torquing mechanism. In this case, disconnect the electrical power to the feed pump, and ensure that the feed tank is properly sealed. Preferably, the tie rod nuts should be loosened and retightened to a torque of 10-15 Nm (88-133 in-lbs). This step must only be done when the system is depressurized and the holder does not contain fluid. It is not recommended to store cassettes in holders that are compressed through a hydraulic sealing mechanism.

Chemical compatibility

Lab and Flow cassettes are resistant to most chemicals typically used in bioprocessing. Because chemical compatibility is dependent on the chemical as well as the exposure time and temperature, specific compatibility questions should be answered experimentally.

Table5. Recommended storage conditions

Duration	Temperature	Chemical
Up to 3 days	Ambient (20-25°C)	Water or 1-3% NaCl
Up to 1 month	Ambient (20-25°C)	0.1 N NaOH
Up to 12 months	Refrigerated (4°C)	0.05 N NaOH

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